



## Aidite A-IS Pro Digital Scanner

Aidite A-IS Pro digital is a fully functional dental 3D-Scanner for professionals. Perfect and precise digitalization for prosthetic cases, orthodontic, and surgical planning.

Small and powerful, it has the fastest speed compared to other dental scanners in this segment.

Perfect for labs and milling centers with high demand, it will help users to expand the production and upgrade productivity.

### Aidite A-IS Pro Digital Scanner Highlights:



Rapid scanning



Full function



Small and lightweight

Parameter	Basic Parameters
Size (L×W×H)	260×270×420mm
Weight	5kg
Resolution	1.3MP
Accuracy	≤0.01mm
Output data format	STL, OBJ
Interface	USB3.0
Power	DC 24V
Light source	Blue light

Scanning time	
Occlusion Scan	8s
Upper and lower jaw Scan	12s
Abutment or Die Scan	3s (less than 8 units)
Impression Scan	58s
Dynamic Scanning	
Scan-Head swings	0-5
Texture Scan	

**Application**  
Abutment teeth, stone models, impressions, triple-tray-impressions, articulator and more.

**Blue light scanning**  
Less affect by external influences.

**Compatible with mainstream dental design software**







